




**EverFlow<sup>®</sup>**

**Paclitaxel Coated PTCA Balloon Catheter**


**Efficient Drug Delivery  
and Proven Efficacy**

*Promise for Heart*



## EverFlow® Design Features

The system comprise of drug delivery specific balloon catheter (VesTraveler®) and polymer-free paclitaxel drug coating.



- ★ **Drug Delivery Specific Balloon Catheter (VesTraveler®):** designed specifically for drug and stent delivery with features of semi-compliant, excellent crossability and precision balloon expansion
- ★ **Special HighBind® Coating Technology:** Uniformly 100% surface coating, precision drug loading and controlled drug release profile
- ★ **Evidence-driven of Medicine:** minimum drug loss during delivery and controlled release profile during balloon expansion
- ★ **Proven Long-last vessel drug retention and efficacy:** effective drug dose remain in vessel up to 30days post balloon expansion.

## EverFlow® Paclitaxel Coated PTCA Balloon Catheter

Drug Delivery Specific Balloon Catheter


Uniformly 100% Balloon Surface Coating

Controlled Drug Release:  
• Better drug Adhesion  
• High Drug Transfer to Vessel







Precision Drug Loading

Proven Long-Last Vessel Drug Retention & Efficacy

**Balloon Expanded**






# Evidence-Driven Device Safety and Efficacy



EverFlow®

SeQuent®




EverFlow®

SeQuent®

**QCA and Histomorphometric comparison between SeQuent & EverFlow in porcine coronary arteries 90days post expansion: there are no significant**

**different.** (LA: lumen Area; IEA: Internal Elastic Area; MA: Media Area; NIA: Neointima Area; RFD: vessel reference diameter; MLD: minimal lumen diameter; Lumen loss(%DS)=LLL/MLD post x 100% (inside Data).

## Evidence–Driven Device Safety and Efficacy



**Coronary Vessel Paclitaxel Residual Dosage:** Sufficient and remain up to 35 days post balloon expansion ( Paclitaxel safe & efficacy concentration: 0.0085–0.85ug/ml. Axel et al: Circulation. 1997; 96:636–645; Sollot et al; J. Clin.Invest. 1995; 95: 1869–1876 )

EverFlow® PTX Pharmacokinetics Study (PK)					
Tissues	Time Points(ug/g)				
	0	1 h	24 h	7 days	35 days
Blood	2.62*	0.0084±0.0043 (n=8)	0	0	0
Vessels (Balloon Injured)	/	6.04±1.78 (n=4)	3.64±2.21 (n=5)	0.024±.0082 (n=3)	0.0043±0.0018 (n=3)
Core Life Organs	Myocardium	0.21	0.006	0	0
	Liver	0	0	0	0
	Spleen	0	0	0	0
	Lung	0	0	0.0091	0
	Kidney	0	0	0	0
	Brain	0	0	0	0

**EverFlow® PK Data:** Paclitaxel (PTX) was detected only within one week in the core organs post balloon expansion, and within safe range (Inside data).

# Order Information

EverFlow® Paclitaxel-Coated PTCA Balloon Catheter (sizes)								
Balloon Diameter (mm)	Balloon Length (mm)							
	10	14	18	22	26	30	34	38
2.00	2.00X10	2.00X14	2.00X18	2.00X22	2.00X26	2.00X30	2.00X34	2.00X38
2.25	2.25X10	2.25X14	2.25X18	2.25X22	2.25X26	2.25X30	2.25X34	2.25X38
2.50	2.50X10	2.50X14	2.50X18	2.50X22	2.50X26	2.50X30	2.50X34	2.50X38
2.75	2.75X10	2.75X14	2.75X18	2.75X22	2.75X26	2.75X30	2.75X34	2.75X38
3.00	3.00X10	3.00X14	3.00X18	3.00X22	3.00X26	3.00X30	3.00X34	3.00X38
3.25	3.25X10	3.25X14	3.25X18	3.25X22	3.25X26	3.25X30	3.25X34	3.25X38
3.50	3.50X10	3.50X14	3.50X18	3.50X22	3.50X26	3.50X30	3.50X34	3.50X38
3.75	3.75X10	3.75X14	3.75X18	3.75X22	3.75X26	3.75X30	3.75X34	3.75X38
4.00	4.00X10	4.00X14	4.00X18	4.00X22	4.00X26	4.00X30	4.00X34	4.00X38

**CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

**EverFlow® Indications:**

- ◆ De novo lesions (primary use in the case of stenosis or occlusion) including small vessels (SVD)
- ◆ Residual stenosis after balloon or stent PTCA
- ◆ Pre and post dilatation during coronary stent implantation
- ◆ Acute or impending vascular occlusion

**Instruction for Use:**

Material required to perform a PTCA balloon catheter procedure:

- ◆ Rapid Exchange PTCA catheter
- ◆ Guiding catheter
- ◆ Inflation syringe with gauge
- ◆ Introduction lock set with dilator
- ◆ 0.014 inch × 175 cm guide wire
- ◆ Contrast medium
- ◆ Y–connector with haemostatic adapter

Prior to use, it is essential that all the instruments, catheters and individual components are carefully checked. An integrity test of the balloon (with the protective balloon sheath left in place), a complete removal of air from the system and a thorough control of the connector seals are mandatory.

Detail description for Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device